


TXNDC5 Antibody (R35557)


Catalog No.	Formulation	Size
R35557-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q8NBS9
Gene ID	81567
Applications	Western Blot : 0.05-0.2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This TXNDC5 antibody is available for research use only.



Western blot testing of human lymph node lysate with TXNDC5 antibody at 0.05ug/ml. Isoforms of ~48kDa and ~36kDa may be detected.



Western blot testing of mouse lymph node (A), stomach (B) and intestine (C) lysates with TXNDC5 antibody at 0.1ug/ml. Predicted molecular weight (mouse): ~46kDa.

Description

Additional name(s) for this target protein: Thioredoxin domain containing 5

Application Notes

Optimal dilution of the TXNDC5 antibody should be determined by the researcher.

Immunogen

Amino acids RDGKKVDQYKGKRD were used as the immunogen for this TXNDC5 antibody. This sequence is common to human isoforms 1 and 2.

Storage

Aliquot and store the TXNDC5 antibody at -20oC.