


## UBC3B Antibody (R33159)

| Catalog No.  | Formulation   | Size   |
|--------------|---|--------|
| R33159-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |  |
|----------------------|--|
| Availability         | 1-3 business days  |
| Species Reactivity   | Human, Rat   |
| Predicted Reactivity | Mouse  |
| Format               | Antigen affinity purified  |
| Clonality            | Polyclonal (goat origin)   |
| Isotype              | Goat Ig  |
| Purity               | Antigen affinity   |
| Gene ID              | 54926  |
| Applications         | Immunohistochemistry (FFPE) : 2-4ug/ml<br>Immunoelectron Microscopy : stains chromatoid bodies in rat spermatids<br>ELISA (peptide) LOD : 1:4000 |
| Limitations          | This UBC3B antibody is available for research use only.  |



IHC testing of FFPE human testis with UBC3B antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

### Description

Additional name(s) for this target protein: Ubiquitin-conjugating enzyme UBC3B; UBE2R2

### Application Notes

Optimal dilution of the UBC3B antibody should be determined by the researcher.

### Immunogen

Amino acids AQQQMTSSQKALM were used as the immunogen for this UBC3B antibody.

## **Storage**

Aliquot and store the UBC3B antibody at -20oC.