


UBR2 Antibody (R36493)

Catalog No.	Formulation	Size
R36493-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q8IWW8
Localization	Nuclear, cytoplasmic
Applications	IHC (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This UBR2 antibody is available for research use only.



IHC testing of FFPE human brain (cortex) tissue with UBR2 antibody at 5ug/ml. HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Alternate name(s) for this target: Ubiquitin protein ligase E3 component n-recognin 2; C6orf133; ubiquitin ligase E3 alpha-II

Application Notes

Optimal dilution of the UBR2 antibody should be determined by the researcher.

Immunogen

Amino acids SCPKSGGDKSRAPT were used as the immunogen for this UBR2 antibody.

Storage

Aliquot and store the UBR2 antibody at -20oC.