


UCHL3 Antibody [clone 26U65] (RQ8920)

Catalog No.	Formulation	Size
RQ8920	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human, Mouse, Rat
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	26U65
Purity	Affinity purified
UniProt	P15374
Localization	Nucleus, Cytoplasm, Golgi
Applications	Western Blot : 1:500 Immunofluorescence : 1:50
Limitations	This UCHL3 antibody is available for research use only.



Immunofluorescent staining of FFPE human U-2 OS cells with UCHL3 antibody (green) and Alpha Tubulin mAb (red). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of 1) human K562, 2) human 293T, 3) human A431, 4) human Caco-2, 5) rat heart and 6) mouse heart tissue lysate with UCHL3 antibody. Predicted molecular weight ~26 kDa.

