


UCP2 Antibody (R33648)


Catalog No.	Formulation	Size
R33648-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7351
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This UCP2 antibody is available for research use only.



Western blot testing of rat adipose lysate with UCP2 antibody at 1ug/ml. Predicted molecular weight ~34 kDa.



Western blot testing of 1) rat adipose and 2) mouse spleen lysate with UCP2 antibody at 1ug/ml. Predicted molecular weight ~34 kDa.

Description

Additional name(s) for this target protein: Uncoupling protein 2

Application Notes

Optimal dilution of the UCP2 antibody should be determined by the researcher.

Immunogen

Amino acids DSVKQFYTKGSEH were used as the immunogen for this UCP2 antibody.

Storage

Aliquot and store the UCP2 antibody at -20oC.