


UNC45B Antibody (R35760)

Catalog No.	Formulation	Size
R35760-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q8IWX7
Gene ID	146862
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This UNC45B antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with UNC45B antibody at 0.3ug/ml. Predicted molecular weight: ~104/95kDa (isoforms 1/2).

Description

Additional name(s) for this target protein: Unc-45 homolog B

Application Notes

Optimal dilution of the UNC45B antibody should be determined by the researcher.

Immunogen

Amino acids RSDKLRQKIFKER were used as the immunogen for this UNC45B antibody.

Storage

Aliquot and store the UNC45B antibody at -20oC.