


## UNC5B Antibody (R35840)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R35840-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                          |
|----------------------|----------------------------------------------------------|
| Availability         | 1-3 business days                                        |
| Species Reactivity   | Human                                                    |
| Predicted Reactivity | Cow, Dog, Mouse, Pig, Rat                                |
| Format               | Antigen affinity purified                                |
| Clonality            | Polyclonal (goat origin)                                 |
| Isotype              | Goat Ig                                                  |
| Purity               | Antigen affinity                                         |
| Gene ID              | 219699                                                   |
| Applications         | Western Blot : 2-4ug/ml<br>ELISA (peptide) LOD : 1:64000 |
| Limitations          | This UNC5B antibody is available for research use only.  |



Western blot testing of human bone marrow lysate with UNC5B antibody at 2ug/ml with [B] and without [A] blocking/immunizing peptide. Predicted molecular weight: 104kDa.

## Description

Additional name(s) for this target protein: Unc-5 homolog B

## Application Notes

Optimal dilution of the UNC5B antibody should be determined by the researcher.

## Immunogen

Amino acids SARDWIFQLKTQAHQ were used as the immunogen for this UNC5B antibody.

## Storage

Aliquot and store the UNC5B antibody at -20oC.