


UNC5C Antibody (R34772)


Catalog No.	Formulation	Size
R34772-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8633
Localization	Cytoplasmic
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml Immunofluorescence : suitable ELISA (peptide) LOD : 1:32000
Limitations	This UNC5C antibody is available for research use only.



IHC/IF testing of FFPE islets in mouse pancreas with UNC5C antibody at 50ug/ml (red) and DAPI nuclear counterstain (blue). HIER: citrate buffer.



IHC staining of FFPE human cerebral cortex with UNC5C antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Unc-5 homolog C

Application Notes

Optimal dilution of the UNC5C antibody should be determined by the researcher.

Immunogen

Amino acids NCTVSEEPTGID were used as the immunogen for this UNC5C antibody.

Storage

Aliquot and store the UNC5C antibody at -20oC.