


USP10 Antibody (RQ7110)

Catalog No.	Formulation	Size
RQ7110	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q14694
Applications	Western Blot : 0.5-1ug/ml Direct ELISA : 0.1-0.5ug/ml
Limitations	This USP10 antibody is available for research use only.



Western blot testing of 1) human HeLa, 2) human 293T, 3) human MCF7, 4) human A549, 5) human K562, 6) human SiHa, 7) rat brain, 8) mouse brain and 9) mouse NIH 3T3 cell lysate with USP10 antibody. Predicted molecular weight ~88 kDa but commonly observed at up to ~110 kDa.

Description

Ubiquitin specific peptidase 10, also known as USP10, is an enzyme which in humans is encoded by the USP10 gene. Ubiquitin is a highly conserved protein that is covalently linked to other proteins to regulate their function and degradation. This gene encodes a member of the ubiquitin-specific protease family of cysteine proteases. The enzyme specifically cleaves ubiquitin from ubiquitin-conjugated protein substrates. The protein is found in the nucleus and cytoplasm. It functions as a co-factor of the DNA-bound androgen receptor complex, and is inhibited by a protein in the Ras-GTPase pathway. The human genome contains several pseudogenes similar to this gene. Several transcript variants, some protein-coding and others not protein-coding, have been found for this gene.

Application Notes

Optimal dilution of the USP10 antibody should be determined by the researcher.

Immunogen

Recombinant human protein (amino acids M1-L798) was used as the immunogen for the USP10 antibody.

Storage

After reconstitution, the USP10 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.