


USP28 Antibody (R36098)

Catalog No.	Formulation	Size
R36098-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57646
Applications	Immunohistochemistry (FFPE) : 2.5-5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This USP28 antibody is available for research use only.



IHC testing of FFPE human testis with USP28 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Ubiquitin specific peptidase 28

Application Notes

Optimal dilution of the USP28 antibody should be determined by the researcher.

Immunogen

Amino acids SVELKHYIQEDN were used as the immunogen for this USP28 antibody.

Storage

Aliquot and store the USP28 antibody at -20oC.