


Uteroglobin Antibody / SCGB1A1 [clone AEDA-19] (RQ5374)

Catalog No.	Formulation	Size
RQ5374	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AEDA-19
Purity	Affinity purified
UniProt	P11684
Applications	Western Blot : 1:500-1:2000 Flow Cytometry : 1:50
Limitations	This Uteroglobin antibody is available for research use only.



Western blot testing of human A549 cell lysate with Uteroglobin antibody. Expected molecular weight: 10-16 kDa.

Description

The SCGB1A1 gene encodes a member of the secretoglobin family of small secreted proteins. Uteroglobin has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma. [RefSeq]

Application Notes

Optimal dilution of the Uteroglobin antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human Uteroglobin / SCGB1A1 was used as the immunogen for the Uteroglobin antibody.

Storage

Store the Uteroglobin antibody at -20oC.