


V5 Tag Antibody [clone DA-22] (RQ4680)

Catalog No.	Formulation	Size
RQ4680	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	V5 Tag
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	DA-22
Purity	Affinity purified
Applications	Western Blot : 1:3000-1:10000
Limitations	This V5 Tag antibody is available for research use only.



Western blot testing of a V5 Tag fusion protein lysate with V5 Tag antibody.

Description

Epitope tags consisting of short sequences recognized by well-characterized antibodies have been widely used in the study of protein expression in various systems. The V5 epitope tag is composed of a 14 residue peptide, GKPIPNPLLGLDST, derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5).

Application Notes

Optimal dilution of the V5 Tag antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to V5 was used as the immunogen for the V5 Tag antibody.

Storage

Store the V5 Tag antibody at -20oC.