


VAP-1 Antibody (R33999)

Catalog No.	Formulation	Size
R33999-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8639
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This VAP-1 antibody is available for research use only.



Western blot testing of A549 lysate with VAP-1 antibody at 0.1ug/ml. Predicted molecular weight ~85 kDa but may be observed at higher molecular weights due to glycosylation.

Description

Additional name(s) for this target protein: Vascular adhesion protein 1; AOC3

Application Notes

Optimal dilution of the VAP-1 antibody should be determined by the researcher.

Immunogen

Amino acids ESQAPKTIRDAFC were used as the immunogen for this VAP-1 antibody.

Storage

Aliquot and store the VAP-1 antibody at -20oC.