


VASP Antibody / Vasodilator stimulated phosphoprotein (R35662)


Catalog No.	Formulation	Size
R35662-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7408
Applications	Western Blot : 0.01-0.1ug/ml (human); 1-2ug/ml (mouse) ELISA (peptide) LOD : 1:128000
Limitations	This VASP antibody is available for research use only.



Western blot testing of human platelet lysate with VASP antibody at 0.01ug/ml. Expected molecular weight: ~40/46-50 kDa (unmodified/phosphorylated).



Western blot testing of mouse thymus lysate with VASP antibody at 1ug/ml. Observed molecular weight: ~46/50kDa (unmodified/phosphorylated).

Description

Additional name(s) for this target protein: Vasodilator-stimulated phosphoprotein

Application Notes

Optimal dilution of the VASP antibody should be determined by the researcher.

Immunogen

Amino acids RKVSKQEEASGGPT were used as the immunogen for this VASP antibody.

Storage

Aliquot and store the VASP antibody at -20oC.