


## VDBP Antibody (R33960)

Catalog No.	Formulation	Size
R33960-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2638
Applications	Western Blot : 0.5-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This VDBP antibody is available for research use only.



Western blot testing of human lung lysate with VDBP antibody at 1ug/ml. Predicted molecular weight: ~53kDa.

## Description

Additional name(s) for this target protein: Group-specific component (vitamin D binding protein); GC

## Application Notes

Optimal dilution of the VDBP antibody should be determined by the researcher.

## Immunogen

Amino acids ERGRDYEKNKVCK were used as the immunogen for this VDBP antibody.

## Storage

Aliquot and store the VDBP antibody at -20oC.