


VDBP Antibody / Vitamin D binding protein [clone AOIA-7] (RQ5013)

Catalog No.	Formulation	Size
RQ5013	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AOIA-7
Purity	Affinity purified
UniProt	P02774
Applications	Western Blot : 1:500-1:2000
Limitations	This VDBP antibody is available for research use only.



Western blot testing of human plasma lysate with VDBP antibody. Predicted molecular weight ~53 kDa.

Description

Vitamin D binding protein is involved in vitamin D transport and storage, scavenging of extracellular G-actin, enhancement of the chemotactic activity of C5 alpha for neutrophils in inflammation and macrophage activation.

Application Notes

Optimal dilution of the VDBP antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human VDBP / GC was used as the immunogen for the VDBP antibody.

Storage

Store the VDBP antibody at -20oC.