


VDR Antibody / Vitamin D receptor (R33239)

Catalog No.	Formulation	Size
R33239-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7421
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This VDR antibody is available for research use only.



Western blot testing of human brain lysate with VDR antibody at 0.3ug/ml. Predicted molecular weight 48/54 kDa (isoforms 1/2).

Description

Additional name(s) for this target protein: Vitamin D receptor

Application Notes

Optimal dilution of the VDR antibody should be determined by the researcher.

Immunogen

Amino acids CGNQDYKYRVSD were used as the immunogen for this VDR antibody.

Storage

Aliquot and store the VDR antibody at -20oC.