


VGF Antibody (R31090)

Catalog No.	Formulation	Size
R31090	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal
UniProt	O15240
Applications	Western Blot : 0.5-1ug/ml
Limitations	This VGF antibody is available for research use only.



Western blot testing of VGF antibody and Lane 1: human U87; 2: (h) SHG-44; 3: mouse Neuro-2a cell lysate. Expected molecular weight: ~90/70kDa (unprocessed/processed).

Description

VGF, also called NERP-2 (Neuroendocrine regulatory peptide-2), is a protein and neuropeptide that may play a role in regulating energy homeostasis, metabolism. Canu et al. (1997) noted that rat Vgf encodes a predicted 70-kD polypeptide that shares similarities with the secretogranin/chromogranin family and is found in the secretory granules of subsets of neurons and endocrine cells. By fluorescence in situ hybridization, Canu et al. (1997) assigned the VGF gene to 7q22. Canu et al. (1997) demonstrated that the single-copy human gene spans 6 kb of genomic DNA and contains 2 exons. The entire protein is encoded by exon 2, while exon 1 contains only 5-prime untranslated sequence. The structural organization of the human gene is similar to that described for the rat gene (Salton et al., 1991), and both the translated and the untranslated regions show a high degree of sequence homology to the rat gene.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the VGF antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

An amino acid sequence from the C-terminus of human VGF (QARRAQEEAEAEERRRLQEQ) was used as the immunogen for this VGF antibody.

Storage

After reconstitution, the VGF antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.