


Vimentin Antibody (R35503)


Catalog No.	Formulation	Size
R35503-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7431
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This Vimentin antibody is available for research use only.



Western blot testing of HeLa [A] and Jurkat [B] lysates with Vimentin antibody at 2ug/ml. Predicted molecular weight ~54 kDa.



Western blot testing of 1) mouse ovary (Ab at 0.1ug/ml) and 2) rat ovary (Ab at 2ug/ml) lysates with Vimentin antibody. Predicted molecular weight ~54 kDa.

Description

Additional name(s) for this target protein: VIM

Application Notes

Optimal dilution of the Vimentin antibody should be determined by the researcher.

Immunogen

Amino acids QVINETSQHHDDLE were used as the immunogen for this Vimentin antibody.

Storage

Aliquot and store the Vimentin antibody at -20oC.