


Vitamin D-binding protein Antibody (R33212)

Catalog No.	Formulation	Size
R33212-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2638
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This Vitamin D-binding protein antibody is available for research use only.



Western blot testing of lung lysate with Vitamin D-binding protein antibody at 0.03ug/ml.
Predicted molecular weight: ~53kDa.

Description

Additional name(s) for this target protein: Group-specific component (vitamin D binding protein); GC

Application Notes

Optimal dilution of the Vitamin D-binding protein antibody should be determined by the researcher.

Immunogen

Amino acids CDNLSTKNSKFED were used as the immunogen for this Vitamin D-binding protein antibody.

Storage

Aliquot and store the Vitamin D-binding protein antibody at -20oC.