


von Willebrand Factor Antibody / VWF [clone IHF-22] (RQ5050)

Catalog No.	Formulation	Size
RQ5050	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	IHF-22
Purity	Affinity purified
UniProt	P04275
Applications	Western Blot : 1:500-1:2000
Limitations	This von Willebrand Factor antibody is available for research use only.



Western blot testing of human serum lysate with von Willebrand Factor antibody.
Predicted molecular weight ~309 kDa.

Description

VWF acts as a chaperone for coagulation factor VIII, delivering it to the site of injury, stabilizing its heterodimeric structure and protecting it from premature clearance from plasma. [UniProt]

Application Notes

Optimal dilution of the von Willebrand Factor antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human Von Willebrand Factor / VWF was used as the immunogen for the von Willebrand

Factor antibody.

Storage

Store the von Willebrand Factor antibody at -20oC.