


VPS16 Antibody (R33366)


Catalog No.	Formulation	Size
R33366-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64601
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This VPS16 antibody is available for research use only.



Western blot testing of human HepG2 lysate with VPS16 antibody at 1ug/ml. Predicted molecular weight: ~95/78kDa (isoforms 1/2).



Western blot testing of mouse spleen lysate with VPS16 antibody at 0.5ug/ml. Predicted molecular weight: ~95/78kDa (isoforms 1/2).

Description

Additional name(s) for this target protein: Vacuolar protein sorting 16

Application Notes

Optimal dilution of the VPS16 antibody should be determined by the researcher.

Immunogen

Amino acids CYKVQQKDVSEDED were used as the immunogen for this VPS16 antibody.

Storage

Aliquot and store the VPS16 antibody at -20oC.