


WNT2 Antibody [clone CAG-23] (RQ5276)

Catalog No.	Formulation	Size
RQ5276	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	CAG-23
Purity	Affinity purified
UniProt	P09544
Applications	Western Blot : 1:500-1:2000
Limitations	This WNT2 antibody is available for research use only.



Western blot testing of 1) human HeLa and 2) human SK-BR-3 cell lysate with WNT2 antibody. Predicted molecular weight ~40 kDa.

Description

The WNT2 gene is a member of the WNT gene family. The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. [RefSeq]

Application Notes

Optimal dilution of the WNT2 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human WNT2 was used as the immunogen for the WNT2 antibody.

Storage

Store the WNT2 antibody at -20oC.