


WNT3 Antibody [clone 4F11-H1-F4] (F54041)

Catalog No.	Formulation	Size
F54041-0.1ML	In PBS with 50% glycerol and 0.03% ProClin 300	0.1 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1
Clone Name	4F11-H1-F4
Purity	Protein G affinity
UniProt	P56703
Applications	Western Blot : 1:1000
Limitations	This WNT3 antibody is available for research use only.



Western blot of Wnt3 A) transfected and (B) untransfected CHO-K1 cell lysate with WNT3 antibody at 1:1000. Predicted molecular weight 39~65 kDa depending on glycosylation level.

Description

Ligand for members of the frizzled family of seven transmembrane receptors. Wnt-3 and Wnt-3a play distinct roles in cell-cell signaling during morphogenesis of the developing neural tube.

Application Notes

The stated application concentrations are suggested starting points. Titration of the WNT3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A human recombinant partial protein was used as the immunogen for this WNT3 antibody.

Storage

Store the WNT3 antibody at -20oC.