

WNT3A Antibody (R31233)

Catalog No.	Formulation	Size
R31233	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal
UniProt	P56704
Applications	Western Blot : 0.5-1ug/ml
Limitations	This WNT3A antibody is available for research use only.



Western blot testing of WNT3A antibody and SMMC-7721 lysate. Expected/observed size 35~44KD depending on glycosylation level.

Description

Wingless-type MMIV integration site family member 3A is a protein that in humans is encoded by the WNT3A gene. This gene is a member of the WNT family. The gene was assigned to human chromosome 1q42.13. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. It encodes a protein showing 96% amino acid identity to mouse Wnt3A protein, and 84% to human Wnt3 protein.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the WNT3A antibody may be required

due to differences in protocols and secondary/substrate sensitivity.

Immunogen

An amino acid sequence from the C-terminus of human Wingless-type MMIV integration site family member 3A (WVETLRPRYTYFKVPT) was used as the immunogen for this WNT3A antibody (100% mouse homology).

Storage

After reconstitution, the WNT3A antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.