


## WNT3A Antibody (R33241)


Catalog No.	Formulation	Size
R33241-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7473
Localization	Cytoplasmic, secreted
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This WNT3A antibody is available for research use only.



Western blot testing of HEK293 lysate with WNT3A antibody at 1ug/ml. Predicted molecular weight 35~44 kDa depending on glycosylation level.



IHC staining of FFPE human breast with WNT3A antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

## Description

Additional name(s) for this target protein: Wingless-type MMTV integration site family, member 3

## Application Notes

Optimal dilution of the WNT3A antibody should be determined by the researcher.

## Immunogen

Amino acids CGRGHNTRTEKRKEK were used as the immunogen for this WNT3A antibody.

## Storage

Aliquot and store the WNT3 antibody at -20oC.