


XPO5 Antibody / Exportin-5 (RQ7434)


Catalog No.	Formulation	Size
RQ7434	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q9HAV4
Localization	Nuclear, cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Immunofluorescence : 5ug/ml Flow Cytometry : 1-3ug/million cells Direct ELISA : 0.1-0.5ug/ml
Limitations	This XPO5 antibody is available for research use only.



Immunofluorescent staining of FFPE human U-2 OS cells with XPO5 antibody (green) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of human 1) Caco-2 and 2) Daudi cell lysate with XPO5 antibody.
Predicted molecular weight ~136 kDa.



Flow cytometry testing of human Daudi cells with XPO5 antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= XPO5 antibody.

Description

Exportin-5 (XPO5) is a protein that in humans is encoded by the XPO5 gene. The International Radiation Hybrid Mapping Consortium mapped the XPO5 gene to chromosome 6. This gene encodes a member of the karyopherin family that is required for the transport of small RNAs and double-stranded RNA-binding proteins from the nucleus to the cytoplasm. The encoded protein translocates cargo through the nuclear pore complex in a RanGTP-dependent process.

Application Notes

Optimal dilution of the XPO5 antibody should be determined by the researcher.

Immunogen

E. coli-derived recombinant human protein (amino acids M1-L1104) was used as the immunogen for the XPO5 antibody.

Storage

After reconstitution, the XPO5 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.