


ZAP70 Antibody (R35850)

Catalog No.	Formulation	Size
R35850-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7535
Applications	Western Blot : 0.3-1ug/ml Immunoprecipitation : suitable ELISA (peptide) LOD : 1:32000
Limitations	This ZAP70 antibody is available for research use only.



Western blot testing of Jurkat lysate with ZAP70 antibody at 0.3ug/ml.
Predicted/observed molecular weight: 70kDa.

Description

Additional name(s) for this target protein: Zeta-chain associated protein kinase 70kDa

Application Notes

Optimal dilution of the ZAP70 antibody should be determined by the researcher.

1. This ZAP70 antibody is specific for isoform 1.

Immunogen

Amino acids SDGYTPEPARITSPD were used as the immunogen for this ZAP70 antibody.

Storage

Aliquot and store the ZAP70 antibody at -20oC.