


ZDHHC8 Antibody (R35424)

Catalog No.	Formulation	Size
R35424-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	29801
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This ZDHHC8 antibody is available for research use only.



Western blot testing of human brain [frontal cortex] lysate with [B] and without [A] immunizing peptide using ZDHHC8 antibody at 0.1ug/ml. Predicted/observed molecular weight: ~81kDa.

Description

Additional name(s) for this target protein: Zinc finger, DHHC-type containing 8

Application Notes

Optimal dilution of the ZDHHC8 antibody should be determined by the researcher.

1. This ZDHHC8 antibody is specific to isoform 2.

Immunogen

Amino acids QRDHPQLKTPPSK were used as the immunogen for this ZDHHC8 antibody.

Storage

Aliquot and store the ZDHHC8 antibody at -20oC.