


Zebrafish Akt3 Antibody (RZ1120)

Catalog No.	Formulation	Size
RZ1120	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	D9IL79
Applications	Western Blot : 0.5-1 ug/ml
Limitations	This Zebrafish Akt3 antibody is available for research use only.



coordination of mitochondrial biogenesis with growth factor-induced increases in cellular energy demands.

Application Notes

Optimal dilution of the Zebrafish Akt3 antibody should be determined by the researcher.

Immunogen

An E.coli-derived zebrafish Akt3 recombinant protein (amino acids M104-L154) was used as the immunogen for the Zebrafish Akt3 antibody.

Storage

After reconstitution, the Zebrafish Akt3 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.