


Zebrafish Ascl1a Antibody / Achaete-scute homolog 1 (RZ1203)

Catalog No.	Formulation	Size
RZ1203	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q90259
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Ascl1a antibody is available for research use only.



Western blot analysis of Ascl1a protein using Zebrafish Ascl1a antibody and zebrafish embryo tissue lysates. Predicted molecular weight ~22 kDa.