


## Zebrafish Cdkn1a Antibody / p21 / Waf1 / Cip1 (RZ1258)

Catalog No.	Formulation	Size
RZ1258	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

**Bulk quote request**

<b>Availability</b>	2-3 weeks
<b>Species Reactivity</b>	Zebrafish
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity chromatography
<b>Buffer</b>	Lyophilized from 1X PBS with 2% Trehalose
<b>UniProt</b>	A0A8M1RFS1
<b>Applications</b>	Western Blot : 0.5-1ug/ml
<b>Limitations</b>	This Zebrafish Cdkn1a antibody is available for research use only.



Western blot analysis of Cdkn1a protein using Zebrafish Cdkn1a antibody and 1) zebrafish head tissue lysates and 2) zebrafish embryo tissue lysate. Predicted molecular weight ~15 kDa.