


Zebrafish Cyp11a1.1 Antibody / Cyp11a1 (RZ1213)

Catalog No.	Formulation	Size
RZ1213	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q7SYJ6
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Cyp11a1.1 antibody is available for research use only.



Western blot analysis of Cyp11a1.1 protein using Zebrafish Cyp11a1.1 antibody and 1) zebrafish head, 2) whole female zebrafish and 3) whole male zebrafish tissue lysate. Predicted molecular weight ~59 kDa.