


Zebrafish Cyp26a1 Antibody / Cytochrome P450 26A1 (RZ1216)

Catalog No.	Formulation	Size
RZ1216	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P79739
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Cyp26a1 antibody is available for research use only.



Western blot analysis of Cyp26a1 protein using Zebrafish Cyp26a1 antibody and 1) whole female zebrafish tissue lysates and 2) whole male zebrafish tissue lysates. Predicted molecular weight ~56 kDa.