


Zebrafish Estrogen receptor 2b Antibody / Esr2b (RZ1221)

Catalog No.	Formulation	Size
RZ1221	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	B8JL3
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Estrogen receptor 2b antibody is available for research use only.



Western blot analysis of Estrogen receptor 2b protein using Zebrafish Estrogen receptor 2b antibody and 1) zebrafish head, 2) whole female zebrafish, 3) whole male zebrafish and 4) zebrafish embryo tissue lysate. Predicted molecular weight ~66 kDa.