


Zebrafish Gfap Antibody (F52537)


Catalog No.	Formulation	Size
F52537-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52537-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q58EE9
Applications	IHC (Paraffin) : 1:25 Western Blot : 1:1000
Limitations	This Zebrafish Gfap antibody is available for research use only.



IHC analysis of FFPE zebrafish muscle section using Gfap antibody; Ab was diluted at 1:25.



Western blot analysis of lysate from zebrafish brain, zebrafish whole tissue lysate (left to right) using Gfap antibody. Ab was diluted at 1:1000 for each lane.

