


Zebrafish Her4 Antibody / Her4.1 / Her4.2 / Her4.3 / Her4.5 (RZ1188)


| Catalog No. | Formulation | Size |
|-------------|---|--------|
| RZ1188 | 0.5mg/ml if reconstituted with 0.2ml sterile DI water | 100 ug |

Bulk quote request

| | |
|---------------------------|--|
| Availability | 2-3 weeks |
| Species Reactivity | Zebrafish |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (rabbit origin) |
| Isotype | Rabbit Ig |
| Purity | Antigen affinity chromatography |
| Buffer | Lyophilized from 1X PBS with 2% Trehalose |
| UniProt | Q90466 |
| Localization | Nuclear |
| Applications | Immunohistochemistry (FFPE) : 2-5ug/ml |
| Limitations | This Zebrafish Her4 antibody is available for research use only. |



IHC staining of zebrafish Her4 protein using Zebrafish Her4 antibody, HRP-labeled secondary and DAB substrate. Her4 was detected in a paraffin-embedded section of zebrafish brain tissue. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



IHC staining of zebrafish Her4 protein using Zebrafish Her4 antibody. Her4 was detected in a paraffin-embedded section of zebrafish brain tissue. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.

