


Zebrafish Kv1.2 Antibody / Kcna2a (RZ1151)

Catalog No.	Formulation	Size
RZ1151	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	E7F8M2
Applications	Western Blot : 0.5-1 ug/ml
Limitations	This Zebrafish Kv1.2 antibody is available for research use only.



Western blot analysis of Kv1.2/Kcna2a protein using Zebrafish Kv1.2 antibody and 1) zebrafish head, 2) whole female zebrafish, 3) whole male zebrafish and 4) zebrafish embryo tissue lysate. Predicted molecular weight ~57 kDa.