


Zebrafish Mcad Antibody / Acadm (RZ1138)

Catalog No.	Formulation	Size
RZ1138	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	A2CG95
Applications	Western Blot : 0.5-1 ug/ml
Limitations	This Zebrafish Mcad antibody is available for research use only.



Western blot analysis of Acadm/Mcad protein using Zebrafish Mcad antibody and 1) zebrafish head, 2) female zebrafish viscera, 3) male zebrafish viscera, 4) whole zebrafish and 5) zebrafish embryo tissue lysate. Predicted molecular weight ~46 kDa.