


## Zebrafish Mcm4 Antibody (RZ1139)

Catalog No.	Formulation	Size
RZ1139	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

**Bulk quote request**

<b>Availability</b>	2-3 weeks
<b>Species Reactivity</b>	Zebrafish
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity chromatography
<b>Buffer</b>	Lyophilized from 1X PBS with 2% Trehalose
<b>UniProt</b>	Q6NZV2
<b>Applications</b>	Western Blot : 0.5-1 ug/ml
<b>Limitations</b>	This Zebrafish Mcm4 antibody is available for research use only.



Western blot analysis of Mcm4 protein using Zebrafish Mcm4 antibody and zebrafish embryo tissue lysate. The predicted molecular weight of MCM4 is ~97 kDa.