


## Zebrafish Myca Antibody / Myc-A / c-Myc-A (RZ1241)

| Catalog No. | Formulation   | Size   |
|-------------|---|--------|
| RZ1241      | 0.5mg/ml if reconstituted with 0.2ml sterile DI water | 100 ug |

**Bulk quote request**

|                           |  |
|---------------------------|--|
| <b>Availability</b>       | 2-3 weeks  |
| <b>Species Reactivity</b> | Zebrafish  |
| <b>Format</b>             | Antigen affinity purified  |
| <b>Clonality</b>          | Polyclonal (rabbit origin)                                       |
| <b>Isotype</b>            | Rabbit Ig  |
| <b>Purity</b>             | Antigen affinity chromatography                                  |
| <b>Buffer</b>             | Lyophilized from 1X PBS with 2% Trehalose                        |
| <b>UniProt</b>            | P52160   |
| <b>Applications</b>       | Western Blot : 0.5-1ug/ml  |
| <b>Limitations</b>        | This Zebrafish Myca antibody is available for research use only. |



Western blot analysis of Myca protein using Zebrafish Myca antibody and 1) zebrafish head, 2) whole female zebrafish, 3) whole male zebrafish and 4) zebrafish embryo tissue lysate. Predicted molecular weight ~48 kDa.