


## Zebrafish N-cadherin Antibody / Ncad / Cdh2 / Cadherin 2 (RZ1246)

Catalog No.	Formulation	Size
RZ1246	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

**Bulk quote request**

<b>Availability</b>	2-3 weeks
<b>Species Reactivity</b>	Zebrafish
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity chromatography
<b>Buffer</b>	Lyophilized from 1X PBS with 2% Trehalose
<b>UniProt</b>	Q90275
<b>Applications</b>	Western Blot : 0.5-1ug/ml
<b>Limitations</b>	This Zebrafish N-cadherin antibody is available for research use only.



Western blot analysis of N-cadherin protein using Zebrafish N-cadherin antibody and 1) zebrafish head tissue lysates and 2) zebrafish embryo tissue lysate. Predicted molecular weight ~97 kDa but may be observed at higher molecular weights due to glycosylation.