


Zebrafish Nkx2.2a Antibody / Nkx2.2 (RZ1250)

Catalog No.	Formulation	Size
RZ1250	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q90481
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Nkx2.2a antibody is available for research use only.



Western blot analysis of Nkx2.2a protein using Zebrafish Nkx2.2a antibody and 1) zebrafish head, 2) whole female zebrafish and 3) whole male zebrafish tissue lysate. Predicted molecular weight ~30 kDa, commonly observed at 36-40 kDa.