


Zebrafish Npas4l Antibody / Cloche / Neuronal PAS domain-containing protein 4-like (RZ1252)

Catalog No.	Formulation	Size
RZ1252	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P0DOC7
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Npas4l antibody is available for research use only.



Western blot analysis of Npas4l protein using Zebrafish Npas4l antibody and 1) zebrafish head, 2) whole female zebrafish, 3) whole male zebrafish and 4) zebrafish embryo tissue lysate. Predicted molecular weight ~69 kDa.