

Zebrafish Otud5 Antibody / Otud5a / Otud5b (RZ1257)

Catalog No.	Formulation	Size
RZ1257	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q08BW0, F1QPN8
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Otud5 antibody is available for research use only.

