

## Zebrafish Patched Antibody (R35480)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R35480-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                                     |
|--------------------|---------------------------------------------------------------------|
| Availability       | 1-3 business days                                                   |
| Species Reactivity | Zebrafish                                                           |
| Format             | Antigen affinity purified                                           |
| Clonality          | Polyclonal (goat origin)                                            |
| Isotype            | Goat Ig                                                             |
| Purity             | Antigen affinity                                                    |
| UniProt            | Q98864                                                              |
| Applications       | Immunofluorescence : suitable (1)<br>ELISA (peptide) LOD : 1:32000  |
| Limitations        | This Zebrafish Patched antibody is available for research use only. |



### Description

Additional name(s) for this target protein: Zebrafish patched 1, ptc1; ptc1

### Application Notes

Optimal dilution of the Zebrafish Patched antibody should be determined by the researcher.

1. Has been used to study tumors in zebrafish (PMID: 19555497)

## Immunogen

Amino acids QTGSKKEPFNYSQ were used as the immunogen for this Zebrafish Patched antibody.

## Storage

Aliquot and store the Zebrafish Patched antibody at -20oC.