


Zebrafish Pdx1 Antibody / Pancreas/duodenum homeobox protein 1 (RZ1267)

Catalog No.	Formulation	Size
RZ1267	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q6DC85
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Pdx1 antibody is available for research use only.



Western blot analysis of Pdx1 protein using Zebrafish Pdx1 antibody and 1) zebrafish head, 2) whole female zebrafish, 3) whole male zebrafish tissue lysate. Predicted molecular weight ~28 kDa but may be observed at higher molecular weights due to post-translational modifications.