


Zebrafish Pgam1 Antibody / Pgam1a / Pgam1b / Phosphoglycerate mutase (RZ1268)

Catalog No.	Formulation	Size
RZ1268	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q7SZR4, B8A4H6
Applications	Immunohistochemistry (FFPE) : 2-5ug/ml
Limitations	This Zebrafish Pgam1 antibody is available for research use only.



IHC staining of zebrafish Pgam1a/b protein using Zebrafish Pgam1 antibody, HRP-labeled secondary and DAB substrate. Pgam1 was detected in a paraffin-embedded section of zebrafish heart tissue. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.