


Zebrafish Pi4kb Antibody / Phosphatidylinositol 4-kinase beta (RZ1270)


Catalog No.	Formulation	Size
RZ1270	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q49GP3
Localization	Cytoplasm (ER, Mitochondria)
Applications	Immunohistochemistry (FFPE) : 2-5ug/ml
Limitations	This Zebrafish Pi4kb antibody is available for research use only.



IHC staining of zebrafish Pi4kb protein using Zebrafish Pi4kb antibody, HRP-labeled secondary and DAB substrate. Pi4kb was detected in a paraffin-embedded section of zebrafish spinal tissue. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



IHC staining of zebrafish Pi4kb protein using Zebrafish Pi4kb antibody, HRP-labeled secondary and DAB substrate. Pi4kb was detected in a paraffin-embedded section of zebrafish brain tissue. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.

