


Zebrafish Polr2a Antibody / DNA-directed RNA polymerase subunit (RZ1272)

Catalog No.	Formulation	Size
RZ1272	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	2-3 weeks
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	A0A0R4IMS9
Applications	Western Blot : 0.5-1ug/ml
Limitations	This Zebrafish Polr2a antibody is available for research use only.



Western blot analysis of Polr2a protein using Zebrafish Polr2a antibody and zebrafish embryo tissue lysates. Predicted molecular weight ~217 kDa but commonly observed at up to ~270 kDa (human similarity).